# Making any gene targetable The **OMNI<sup>TM</sup>** panel of novel engineered nucleases **answers key needs** in genome editing

Liat Rockah, Ella Segal, Bar Ben Baruch, Rafi Emmanuel, Lior Izhar

## The **OMNI<sup>TM</sup> nuclease** panelmaking any gene targetable

Successful gene editing requires developing a carefully optimized. Using this powerful platform, we have generated a composition per target gene and target organ. Access to a variety of panel of novel nucleases OMNI<sup>TM</sup> that are used in nucleases and guides, characterized by high activity and specificity and our own clinical programs, unlocking the full applicable PAM, coupled with the right delivery tool is key. Different potential of genome editing. These Type II nucleases editing objectives require specific nuclease characteristics, such as are diverse in size and compatible with all delivery nuclease size or PAM usage, while exhibiting superb specificity. Emendo modalities; PAM usage diversity allows for 86% Biotherapeutics has developed a dual platform technology combining a genome coverage; and highly specific engineered discovery pipeline and cutting-edge protein-engineering capabilities, variants enable allele specific editing. Together all of supported by extensive computational and machine learning tools.

the above makes any gene targetable.



## Compatible with all delivery modalities

### Short AAV-delivery-compatible nucleases - OMNI-A2 (~1,050aa)

nucleases

Adeno Associate Virus (AAV) based vectors are an important part in the delivery toolbox but have limited payload capacity. Using its discovery platform, Emendo Biotherapeutics has identified short CRISPR nucleases compatible with AAV packaging. Protein engineering was used to further improve the editing efficiency of these short nucleases, as exemplified by OMNI-A2 and its hyper active variant OMNI-A2-V20.





OMNI-A2 and engineered OMNI-A2-V20 Highly active, short, AAV packaging compatible novel



AAV delivered OMNI-A2 editing



## Safety

OMNI-A1 was identified using Emendo Biotherapeutics' nuclease discovery platform, exhibiting high editing activity as well as low off target effects. Due to the key safety aspect of off-target effects, OMNI-A1 was optimized using the engineering platform to eliminate off-target editing. CIRCLEseq analysis found OMNI-A1-V10, an engineered variant of OMNI-A1, to have significantly less off-target effects than the wild-type nuclease or spCas9. Importantly, no Off-Target Mediated Translocations (OMTs) were detected using unbiased CASTseq analysis of OMNI-A1-V10 in target cells.

## Genome accessibility

### **OMNI<sup>TM</sup>** Nuclease panel to overcome PAM Constraints

The diversity of OMNI<sup>TM</sup> nucleases particularly PAM site diversity widens significantly genome accessibility and enables targeting of genomic sites that are not accessible using NGG nucleases. Accumulatively the OMNI<sup>TM</sup> panel covers approximately 86% of the genome, making any gene targetable.

Short nucleases					
OMNI™	PAM	Length (aa)			
OMNI-A2	NGGNNNNN	1062			
OMNI-A7	NNGNRCNN	1097			
OMNI-A8	NYAGCNNN	1109			
OMNI-A13	NNNNCMAN	1091			
OMNI-A31	NNGRVNNN	1054			
OMNI-A34	NNNNCVTA	1078			

#### Long nucleases

OMNI™	PAM	Length (aa)
OMNI-A1	NGGNNNNN	1370
OMNI-A3	NNGNRANN	1154
OMNI-A54	NRRNATNN	1107
OMNI-A4	NNRACTTN	1348
OMNI-A8	NYAGCNNN	1109
OMNI-A21	NRGGGCRN	1215

#### Major IP issue in the cell therapy space: heavily patented NGG guides

A little recognized intellectual property (IP) limitation is that of patented guide sequences. For many gene targets, for example CAR-T gene targets, multiple patents by various entities cover applicable NGG guide sequences, making them unavailable for companies that desire using CRISPR as their editing tool. Using OMNI non-NGG nucleases and guides overcomes guide IP barriers.

Gene name	<pre># of patent families directed to     guide sequences</pre>	<pre># of NGG guides left after patent screen</pre>	# of guides for non-NGG nucleases selected for Emendo Biotherapeutics' screening
TRAC	56	0/20	15
TRBC1	23	1/20	3
TRBC2	23	1/21	5
CD3e	4	10/79	18
B2M	42	O/33	12
CIITA	9	0/294	35
PD1	51	0/196	29

### Super-specific and highly active **OMNI<sup>TM</sup>** nucleases

- No off-target effects
- No translocations
- Allele specific editing







MaxEdit 0.98

> 0.81 0.85 0.93 0.88 0.89



### surface proteins



#### CAR-T ready **OMNI<sup>TM</sup>** and guide compositions

Target gene	% Editing ( by NG
B2M	91
CD3	70
CIITA	82
CTLA4	75
FASLG	68
HAVCR2	94
HLA-E	94
PDI	44
TET	88
TRAC	90
CISH	70



Gene knockout using OMNI-A4 results in a dramatic reduction in expression of cell

We are open to collaboration on the **OMNI<sup>TM</sup>** nuclease panel and editing capabilities

For further information please contact u

